




## Quartz Cuvettes With PTFE Caps

The utilization of [Quartz cuvettes with PTFE caps](#) isn't generally vital and it is subject to the frequency range been focused on. The transmission rate for quartz cuvette is high (83% or higher) UV length range while plastic (PS or PMMA) or Optical glass don't. At the point when natural solvents are utilized, glass or quartz is ideal on the grounds that the corvettes show higher opposition when contrasted with others produced using plastic.



# Quartz Cuvettes With PTFE Caps

**Quartz Cuvettes With PTFE Caps provide high-precision optical clarity for spectroscopy, ensuring chemical resistance and secure sealing for accurate laboratory analysis.**



## Contact Us!



+1 877-419-0711



[info@alphananotechne.com](mailto:info@alphananotechne.com)



[www.alphananotechne.com](http://www.alphananotechne.com)