

## **Tubigrip Size Chart Pdf**

File name: Tubigrip Size Chart Pdf
Rating: 4.5/5 (Based on 8367 votes)
29505 downloads
→ Tubigrip Size Chart Pdf

Over 1 Million Products · Super Low PricesA+ Rated – Better Business Bureau (BBB). Directions for use Warning: Contains natural rubber latex. Cut Tubigrip to twice length required for limb, allowing an extra 2–3 cm for overlap. Pull Tubigrip onto limb like a stocking. Double . Top Lymphedema Brands · Top Brands & Products · Free Shipping · Highly AbsorbentBrands: Juzo, Medi, Jobst, BiaCare, Solaris, CircAid, Sigvaris, BSN Medical. Tubigrip® provides tissue support in the treatment of strains and sprains, soft tissue injuries, general edema, post-burn scarring and ribcage injuries. Cut Tubigrip® 10m to twice length required for limb, allowing an extra 2–3cm for overlap. Pull Tubigrip® onto limb like a stocking. Tubigrip is an interwoven stockinette with covered latex rubber yarns providing compression and support used for joint and venous conditions. Low, low-medium, medium and high compression can be obtained, depending on the size used and the size of the limb that it is applied to. Tubigrip sizes range from size A (small wrist) to L (abdomen) and are available in a variety of lengths. It is important to choose the right size for you to get the right amount of compression and support.