



# Gram Stainer

The Automated Gram Stainer provides fast easy and stable staining process for staining clinical specimens according to GRAM staining method by automated system control from fixation to drying. The Automated Gram Stainer features advanced functions and performances over manual staining with simplicity and convenience. It is uniquely designed to provide accurate and precision stain for up to 10 slides (20 slides by adding Optional Tray) loaded in the gram stainer chamber. The atomizer spray used one time prepared fresh reagent to spray on a sample containing slides fixed in a carousel with a different spray nozzle provide precise and accurate throughput.

[https://www.abimed.com/gram-stainer?utm\\_source=justpaste&utm\\_medium=CDA-Selina\\_Carter](https://www.abimed.com/gram-stainer?utm_source=justpaste&utm_medium=CDA-Selina_Carter)

